



# The Bone Grafts Catalogue

Bio-PRF Bone Grafts: Trusted Support for Predictable Bone Regeneration.



## CONTENTS



### Extract-Oss

PAGE 3

A bone grafting material designed specifically for extraction site management. Extract-Oss is a 50/50 mix of cortico-cancellous bone (FDDBA)



### Perio-Oss

PAGE 4

A bone grafting material designed for periodontal regeneration specifically for the treatment of intrabony and furcation defects. This graft is composed entirely of 100% pure demineralized freeze-dried bone allograft (DFDBA).




## EXTRACT-OSS

- A more potent version can be obtained from a 75/25 mineralized/demineralized formulation. The 25% demineralized freeze-dried bone allograft (DFDBA) allows for faster release of BMPs from within the graft while maintaining dimensional ridge stability owing to its 75% composition of mineralized freeze-dried bone allograft (FDBA).

### AVAILABLE SIZES

### PRICING

Extract-Oss-0.5CC	\$64, 10 or more \$59
Extract-Oss-1CC	\$84, 10 or more \$79
Extract-Oss-2CC	\$149, 10 or more \$139
Extract-Oss-MinDemin-0.5CC	\$67, 10 or more \$62
Extract-Oss-MinDemin-1CC	\$88, 10 or more \$83
Extract-Oss-MinDemin-2CC	\$154, 10 or more \$144



## PERIO-OSS

- This graft is composed entirely of 100% pure demineralized freeze-dried bone allograft (DFDBA). This favors faster regeneration of the periodontal ligament, with cementum and bone regeneration attached via Sharpey's fibers.

### AVAILABLE SIZES

### PRICING

Perio-Oss-0.5CC

\$69, 10 or more \$64

Perio-Oss-1CC

\$88, 10 or more \$84

# BIOePRF

Bio-PRF Bone Grafts: Trusted Support for Predictable Bone Regeneration.

SCAN QR

